

Assessing the Clinical Practicability of Artificial Intelligence Assisted Reporting Using the Paris System for Urinary Cytology

PRESENTED BY

Tien-Jen Liu, M.D.



USCAP 113TH ANNUAL MEETING

**BRINGING
EDUCATION
TO LIFE**





Disclosure

I serve as the Co-Founder and Chief Operating Officer at AlxMed, Inc.



Introduction



Introduction

- Despite adherence to guidelines of The Paris System (TPS), variability in microscopic interpretation persists, demanding objective diagnostic methods.
- AIxURO, a deep-learning tool, enhances digital urine cytology by accurately identifying atypical cells in whole-slide images (WSIs), demonstrating clinical validity for bladder cancer diagnosis.
- Incorporating AIxURO into the clinical urine cytology workflow is expected:
 - To streamline diagnostic efficiency and accuracy
 - To reduce subjectivity in TPS categorization
 - To surpass the capabilities of conventional microscopy

Materials and Methods

A panel of experts categorized 116 urine cytology slides according to TPS guidelines:
 86 NHGUC
 12 AUC
 11 SHGUC
 7 HGUC
(Ground Truth)

The Leica Aperio AT2 scanner

Reader:

- Cytopathologist
- Cytologist A
- Cytologist B



**Arm 1
Microscopy**

Diagnosis results

- TPS category
- Time spent

2 weeks
washout



**Arm 2
WSI Review**

Diagnosis results

- TPS category
- Time spent

2 weeks
washout



**Arm 3
AI-Assisted
(AIxURO) Review**

Diagnosis results


- TPS category
- Time spent

Evaluation metrics


- Performances for each TPS category: Sensitivity/Specificity
- Time spent on diagnostics for each TPS category

Materials and Methods

Arm 2: WSI review



Arm 3: AI-assisted (AIxURO) review



Results-Diagnosis

Ground Truth TPS Category		NHGUC (N= 86)			AUC (N=12)			SHGUC (N=11)			HGUC (N=7)		
Reader	Method	Arm 1	Arm 2	Arm 3	Arm 1	Arm 2	Arm 3	Arm 1	Arm 2	Arm 3	Arm 1	Arm 2	Arm 3
Cytopathologist	Sensitivity	96.5%	93.0%	93.0% ↓	25.0%	41.7%	75.0% ↑	54.5%	27.3%	36.4% ↓	100.0%	100.0%	100.0%
	Specificity	76.7%	76.7%	90.0% ↑	94.2%	91.3%	91.3% ↓	98.1%	100.0%	100.0% ↑	98.2%	95.4%	96.3% ↓
Cytologist A	Sensitivity	70.9%	80.2%	60.5% ↓	25.0%	25.0%	83.3% ↑	9.1%	9.1%	54.5% ↑	100.0%	85.7%	100.0%
	Specificity	93.3%	83.3%	100.0% ↑	88.5%	89.4%	74.0% ↓	93.3%	93.3%	96.2% ↑	78.9%	87.2%	90.8% ↑
Cytologist B	Sensitivity	94.2%	84.9%	86.0% ↓	25.0%	25.0%	33.3% ↑	18.2%	9.1%	9.1% ↓	100.0%	57.1%	85.7% ↓
	Specificity	76.7%	73.3%	80.0% ↑	90.4%	80.8%	85.6% ↓	99.0%	96.2%	97.1% ↓	95.4%	97.2%	93.6% ↓


- NHGUC Diagnostics: AIxURO (Arm 3) demonstrated lower sensitivity yet higher specificity compared to Microscopy (Arm 1).
- AUC Diagnostics: Sensitivity increased with AIxURO, albeit with a trade-off in specificity.
- SHGUC Diagnostics: Varied sensitivity among evaluators with AIxURO; Cytopathologist and Cytologist B noted a decrease, whereas Cytologist A showed an increase.
- HGUC Diagnostics: Cytopathologist and Cytologist A maintained sensitivity with AIxURO, but Cytologist B experienced a decrease.

Results- Time

Cytopathologist

Cytologist A

Cytologist B



- The Cytopathologist recorded the longest diagnostic time for AUC and the shortest for HGUC.
- Cytotechnologists required the most time for SHGUC diagnosis, with the least time spent on NHGUC.
- Compared to Microscopy (Arm 1), AIxURO (Arm 3) demonstrated a reduction in diagnostic time across all categories—a time savings over microscopy of 32-45% for the Cytopathologist, 10-50% for Cytologist A, and 31-62% for Cytologist B.

Discussion & Summary

- AIxURO outperformed Microscopy in diagnostic accuracy for AUC, maintaining comparable accuracy across other TPS categories.
- AIxURO significantly reduced the diagnostic time for each TPS category.
- Observer variations in diagnostic accuracy and time were evident among the three methods tested.
- Reviewing WSI without AI assistance did not improve diagnostic accuracy or efficiency.
- The findings highlight the crucial role of AI in improving the clinical application of the Paris System in urine cytology reporting.



Thank You

BRINGING EDUCATION TO LIFE



USCAP 113TH ANNUAL MEETING



**#IAMUSCAP
#USCAP2024**

Evaluating the Feasibility of a Two-Step Workflow in Artificial Intelligence-Assisted Digital Urine Cytology

Tien-Jen Liu¹, Wen-Chi Yang², Shin-Min Huang², Hui Wen Ho³, Wei-Lei Yang¹, Hsing-Ju Wu³, Cheng-Hung Yeh¹, Shih-Wen Hsu¹, Ming-Yu Lin¹ and Pei-Yi Chu^{3*}

¹AlxMed, Inc., Santa Clara, CA, USA; ²Department of Pathology, Show Chwan Memorial Hospital, Changhua, Taiwan; ³Research Assistant Center, Show Chwan Memorial Hospital, Changhua, Taiwan;


*Correspondence to PY Chu ✉ chu.peiyi@msa.hinet.net

Disclosure: Tien-Jen Liu, Wei-Lei Yang, Cheng-Hung Yeh, Shih-Wen Hsu, and Ming-Yu Lin are employees of AlxMed, Inc.

Introduction

- Routine urine cytology screenings for bladder cancer contribute to a significant clinical workload, primarily due to the low frequency of abnormalities, with the majority of cases being negative.
- Cervical cancer screening employs a two-step Pap test workflow to enhance efficiency: cytotechnologists initially assess slides for interpretations, and slides deemed "abnormal" are referred to cytopathologists for further evaluation, thereby optimizing the workflow and reducing the burden on cytopathologists.
- The efficacy of artificial intelligence (AI)-based digital urine cytology tools in improving the performance and efficiency of bladder cancer screenings has been established.
- This study investigates the potential of incorporating a two-step urine cytology workflow, augmented by AI-assisted tools, into bladder cancer screening protocols to optimize case management and diagnostic accuracy (Figure 1).

Figure 1. Two-Step Workflow for Urine Cytology in Cancer Screening



Materials and Methods

- An expert panel diagnosed 116 urine cytology slides, identifying 30 as positive and 86 as negative for bladder cancer. These slides, forming the "ground truth," were digitized into whole-slide images (WSIs).
- The slides/images were evaluated by three readers (two cytotechnologists and one cytopathologist) using three distinct methods: microscopy (Arm 1), WSI review (Arm 2), and AI-assisted review using AlxURO (Arm 3).
- The evaluation initiated with Arm 1, adopting the two-step workflow for bladder cancer screening. Following a two-week washout period, the readers progressed to the subsequent methods.
- For each method, performance metrics such as sensitivity, specificity, positive predictive value (PPV), negative predictive value (NPV), and total diagnostic time were calculated for each reader.
- The outcomes were analyzed to assess the feasibility of adopting the workflow in cancer screening (Table 1).

Results

Table 1. Performance and Interpretation Time Across Three Methods in the Two-Step Workflow

Method	Order of Readers	The Two-step Workflow				Step 1: Cytotechnologist			Step 2: Cytopathologist		Workflow Spent Time (min)	
		Performance of Binary Diagnosis (Positive/Negative)				Interpretation of Total Cases		Total Spent Time (min)	Interpretation of Positive Cases			
		Sensitivity	Specificity	PPV	NPV	Positive	Negative		Positive	Negative		
Arm 1 Microscopy	Cytotech A + Cytopath	76.7%	96.5%	88.5%	92.2%	53	63	130.2	26	27	80.0	210.2
	Cytotech B + Cytopath	76.7%	97.7%	92.0%	92.3%	28	88	197.3	25	3	47.4	244.7
Arm 2 WSI review	Cytotech A + Cytopath	70.0%	92.4%	80.8%	90.0%	42	74	170.3	26	16	56.8	227.1
	Cytotech B + Cytopath	63.3%	95.4%	82.6%	88.2%	35	81	187.9	23	12	55.9	243.8
Arm 3 AI-assisted review (AlxURO)	Cytotech A + Cytopath	90.0%	93.0%	81.8%	96.4%	64	52	29.0	33	31	43.2	72.2
	Cytotech B + Cytopath	76.7%	93.0%	79.3%	92.0%	36	80	87.0	29	7	23.4	110.4

Cytotech: Cytotechnologist; Cytopath: Cytopathologist; PPV: Positive Predictive Value; NPV: Negative Predictive Value

- In Arm 3, Cytotechnologists A and B identified more positive cases (64 and 36, respectively) than in Arm 1 (53 and 28) and Arm 2 (42 and 35).
- Arm 3 enhanced the sensitivity (90.0% vs. 76.7% in Arm 1) and NPV (96.4% vs. 92.2% in Arm 1), highlighting the efficiency of the AI-assisted two-step workflow. The performance of Arm 2 paralleled that of Arm 1.
- The adoption of the two-step workflow in Arm 1 reduced the total diagnostic time to 210.2 (74.4%) and 244.7 (70.0%) minutes, a reduction from the 282.5 and 349.6 minutes incurred by combining Cytotechnologists A or B's evaluations with the Cytopathologist's.
- Notably, the adoption of the two-step workflow in Arm 3 further condensed the diagnostic time to 72.2 and 110.4 minutes, demonstrating marked time savings compared to both Arm 1 and Arm 2.

Conclusion

- AI-assisted review (Arm 3) matched the diagnostic accuracy of both microscopy (Arm 1) and WSI review (Arm 2) while significantly reducing the time required for cancer screening.
- Arm 3 improved cytotechnologists' efficiency by identifying more positive cases, highlighting its potential to enhance the preliminary screening process.
- The lack of performance or time improvements with Arm 2 emphasizes the crucial role of AI in optimizing the screening workflow.
- These findings underline AI's pivotal contribution to digital urine cytology and its potential for integrating a two-step workflow in bladder cancer screening.